---

# Application of AlazarTech PCI/PCIe Digitizer on Optical Coherence Tomography (OCT)

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Email: goldfire@139.com

QQ: 841655201

Scott Engineering Limited  
309, FuHua Building,  
China Academy of Science and Technology Development,  
009, First High-Tech South St., NanShan District  
ShenZhen, GuangDong, PR China

---

## Application Note: 应用案例:

### Contents 目录

Application Note:.....	1
Introduction .....	2
Optical Coherence Tomography .....	3
Fourier Domain OCT .....	4
Swept Source OCT.....	6
High-Speed Waveform Digitizer.....	8
References .....	12
应用案例 .....	1
引言 .....	2
光学相干断层扫描 (OCT) .....	3
傅立叶域 OCT .....	4
扫频光源 OCT .....	6
高速数据采集卡 .....	8
参考书目 .....	12

### Introduction 引言

In the last decade, Optical Coherence Tomography has emerged as an important diagnostic and research tool in several medical fields, including ophthalmology, cardiology, and neurology. OCT is unique as a non-contact, non-invasive imaging modality that can acquire sub-surface high-resolution cross sectional images of biological tissue [1].

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Wavelength Swept Source (SS) OCT implementations dominate in longer wavelength regions, particularly at 1310nm where there are an abundance of fiber optic components from the telecommunications industry. Light at this wavelength penetrates deeper in tissue than shorter wavelengths, but suffers from high water absorption. Applications of high speed OCT imaging in this wavelength range include: ocular anterior chamber, early cancer detection, cardiovascular imaging, and visualization of early developmental biology [2].

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In this application note, the basics of OCT imaging are described from the perspective of wavelength swept imaging. A typical implementation of an externally clocked SS OCT system is introduced, and schematically illustrated. The features on AlazarTech waveform digitizers specifically suited for SS OCT signal acquisition are described, enabling high speed, high performance imaging.

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## Optical Coherence Tomography 光学相干断层扫描

Medical imaging with OCT is described as an optical analog to Ultrasound. Axial cross sectional data (called A-scans) are generated by the interference of light backscattered from a sample with light reflected from a reference. A sequence of A-scans acquired while transversely scanning the beam is called a Brightness (B)-scan. Volumetric images can be constructed by sequentially acquiring cross-sectional images generated through raster scanning, usually using galvanometer scanning mirrors. A simplified version of an OCT system using a fiber-based Michelson interferometer-type configuration is illustrated in Figure 1.

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
Figure 1: Simplified OCT using a fiber-based Michelson interferometer type configuration

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Reference Arm: 参考臂

Sample Arm: 样品臂

Light backscattered from the sample and reference arms are recombined in the fiber coupler and interfere to produce fringes related to the optical path length difference. Using a low-coherence light source with a spectrum,  $S(k)$ , the light returned from the reference and sample arms,  $I_R$  and  $I_S$ , with a path length mismatch of  $\Delta z$ , the signal acquired at the detector can be represented by

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$$I_{det}(k) = S(k)(I_R + I_S + 2\sqrt{I_R I_S} \cos(2\Delta z k)) \quad (1)$$

The source spectrum,  $S(k)$ , is related by the Fourier transform to  $S(\Delta z)$ , which represents the Point Spread Function (PSF) or axial resolution of the system. Assuming a Gaussian profile having a full width half maximum (FWHM) spectral bandwidth  $\Delta\lambda$ , the width of the PSF is given by [3] [4]

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$$l_c = \frac{2(\ln 2) \lambda^2}{\pi \Delta\lambda} \quad (2)$$

where  $\lambda$  is the central wavelength and  $\Delta\lambda$  is the source bandwidth in wavelength units. In this form, the resolution can be seen to depend not only on the bandwidth of the source, but also on the central wavelength. This implicitly indicates that in OCT, the axial resolution is de-coupled from the transverse resolution, which is dependent on the sample arm imaging optics.

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## Fourier Domain OCT 傅立叶域 OCT

Fourier domain (FD) detection in OCT has gained much attention due to the increased sensitivity and acquisition speed compared to the original method, time domain (TD) OCT. FD OCT uses a stationary reference arm reflector, and the optical path length difference of scatterers in the sample arm relative to the position of the reference arm is encoded by the frequency of the interferometric fringes as a function of the source spectrum,  $S(k)$ . In FD OCT, acquiring the optical spectrum of back-reflected light at permits all depths to be acquired simultaneously rather than sequentially. The signal acquired at the detector is a function of the wavenumber given by Equation (1), simplified for the case of a single reflector in the sample. The location of the reflector in the sample can be extracted by taking the Fourier transform of Equation (1) to obtain

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where the first term on the right hand side is the Fourier transform of the source. The first term in the parenthesis is the DC term, and corresponds to the DC spectrum, and the second term in the parenthesis represents the location of the reflector and its complex conjugate. Complex conjugate artifacts arise because, generally, OCT systems do not acquire the phase of the interferometric signal. Various techniques have been investigated to resolve this problem during imaging, effectively doubling the imaging depth of the system.

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For OCT systems, the maximum sensitivity is determined by the shot noise limit of the system. Typically, the reference arm reflectivity is much higher than the sample arm, where excess photon noise will dominate. However, largely attenuating the references arm will cause receiver noise to dominate. Accordingly, to achieve maximum sensitivity, the reference arm power has to be adjusted to a reflectivity in between the two cases. Consequently, the maximum sensitivity is a function of the reflectivity of the sample  $R_s$ , power incident on the sample  $P_s$  and electronic charge  $e$  [5]. For SD OCT, the Signal to Noise Ratio (SNR) is given by [3]

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$$SNR_{sdoct} = \frac{R_s P_s \Delta t}{2e}, \quad (4)$$

where  $\Delta t$  is the line integration time. Equation (5) shows that the SNR is not dependent on the detection bandwidth (as was the case with Time Domain OCT) but proportional to integration time.

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The sensitivity advantages obtained through spectrally resolved detection, state-of-the-art FD OCT systems permit operation at line rates in the hundreds of kilohertz with image quality superior to older, slower systems. Two types of FD OCT are widely used in research, and are differentiated based on the combination of the source and detector used. Spectral Domain (SD) OCT (also called S OCT) uses a low coherence broadband light source and a spectrometer which uses a diffraction grating to disperse the interferometric signal spectrum across a linear array detector. Swept Source (SS) OCT uses a wavelength swept laser source and photodiode balanced receiver as the detector. Whereas SD OCT spatially encodes spectrally resolved interference signal, SS OCT time-encodes the spectral signal as the narrow laser linewidth sweeps through a broad spectrum [6]. Although the detected signal is acquired in a different manner, the interpretation is the same, i.e. the frequency of the interferometric fringes

corresponds to the location of the sample scatterer. The system level advantages of SS OCT over SD OCT are evidenced through the differences in signal acquisition. SS OCT systems can be designed with faster line rates and longer imaging depths than possible with the spectrometer based systems.

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## Swept Source OCT 扫频光源OCT

There many configurations for implementing a wavelength swept laser for SS OCT, but the principle components are a gain medium and a tunable wavelength filter designed in a feedback configuration as shown in Figure 2.

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
Figure 2: Conceptual wavelength swept laser configuration

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
Gain Medium: 增益介质

Output: 输出

A typical SS OCT system (Figure 3) consists of a wavelength swept laser, an optical circulator, a balanced photodiode detector and a high-speed waveform digitizer configured into an interferometric topology.

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**Figure 3: Swept source interferometric topology**

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Reference arm: 参考臂

Sample arm: 样品臂


Light from the laser goes through the optical circulator, and is divided into the sample and reference arm by the fiber coupler. The backreflected light in the reference and sample arm is recombined by the coupler. The detector arm of the interferometer is directed to one input of a balanced photodiode detector. The light travelling backward through the source arm is recovered by the circulator and delivered to the other input of the balanced photodetector [5]. Balanced detection provides better sensitivity over a single detector because it rejects common mode noise and is more efficient in light usage. The resulting interferometric signal has an ac-coupled appearance.

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The swept source emits one wavelength at a time as the filter is tuned through the optical bandwidth of the gain medium. Due to the nature of the various tunable filters, the sweep rate is not necessarily linear in time, and additional steps are required to sample the spectrum  $S(k)$  at even wavenumber intervals. Some wavelength swept lasers provide a digital pixel clock output synchronized to the filter, for those that do not, a common optoelectronic method to

calibrate wavenumber sampling is illustrated in Figure 4. A pixel clock signal which is fed into the waveform digitizer can be generated by using a separate fixed path interferometer, represented by MZI in Figure 3. As the wavelength swept source is tuned through the source spectrum, the fixed interferometer will cross the DC threshold at points evenly spaced in wavenumber. Using this signal as an external pixel clock for the waveform digitizer permits the SS OCT interferogram to be sampled evenly in wavenumber space [6].

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**Figure 4: A fixed path Mach Zehnder interferometer (MZI) is used to create a sinusoidal signal whose zero crossings are evenly spaced in wavenumber. This signal can be used as a pixel clock for digitization of the interferometric signal**

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
## High-Speed Waveform Digitizer 高速数据采集卡

AlazarTech waveform digitizers enable high-speed acquisition and processing of SS OCT data. The ability of AlazarTech digitizers to guarantee no missed trigger events is crucial particularly in the acquisition of Doppler OCT, which requires calculation of phase difference of adjacent A-scans. High-speed data streaming is achieved through the proprietary Auto Direct Memory Access (AutoDMA) circuitry together with the fully asynchronous driver and dual-DMA engine in AlazarTech digitizers. The combination of these features allows tasks to be performed in parallel fashion. The AutoDMA circuitry uses hardware to re-arm and initiate data transfer which removes software latencies and the dual-port memory acts as a very deep FIFO and dual-DMA engine permitting maximum transfer rates to the PC (limited by the motherboard). The AutoDMA and dual-port memory process is illustrated in Figure 5 and Figure 6, respectively.

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
TRIGGERED DATA ACQUISITION USING DUAL-PORT MEMORY

NOTE 1: No Trigger Events Are Missed - Guaranteed  
 NOTE 2: Over 95% of CPU cycles are available for data processing

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Figure 5: Dual-port memory data acquisition process schematic.

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
Figure 6: AutoDMA acquire and transfer process.

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
AlazarTech digitizers provide several flavours of AutoDMA—traditional, No Pre-Trigger

(NPT), continuous and triggered streaming. For OCT applications, traditional and NPT are commonly used with the latter being the recommended method. The traditional AutoDMA acquires both pre-trigger and post-trigger data while NPT AutoDMA only captures post-trigger data in each buffer. NPT and traditional AutoDMA can easily acquire and transfer data to the PC at a maximum sustained transfer rate of the motherboard without causing an overflow. One advantage of NPT over traditional AutoDMA is its ability to use the entire on-board memory before an overflow flag is asserted; traditional AutoDMA is only limited to 512 buffers on-board memory if no PC data transfer has been performed.

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The ports and architecture schematic of an AlazarTech PCI bus digitizer are shown in Figure 7. Ports consist of an optional external clock input (SMA), two analog input channels (A and B), external trigger analog input and an auxiliary I/O.

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
Figure 7: (Left) Digitizer ports. (Right) Basic architecture schematic of a PCI bus digitizer.

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Wavelength swept lasers have several different designs and the output can be a unidirectional sweep or a bidirectional (forward and backward) sweep. Other outputs from a wavelength swept laser are line trigger signals. The line trigger can be generated from an optical reference such as a Fibre Bragg Grating (FBG) or the square waveform triggers can be based on the drive waveform. For swept laser designs with a bidirectional sweep, the user can choose to use the forward sweep, backward sweep or both sweeps by using a combination of the two triggers available. Each AlazarTech board has two trigger engines, wherein channel A or B can be used as an additional trigger source to the existing external trigger. The configurability of the AlazarTech board triggers provides flexibility for working with the different types of wavelength swept lasers.

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The optional External Clock input provides a solution for k-linear SS OCT acquisition. Sampling the OCT signal linearly in k-space using an external pixel clock reduces the overhead of software resampling. The External Clock input on the AlazarTech waveform digitizers can be used with either a digital pixel clock synchronized with the laser sweep (provided by some commercial wavelength swept lasers) or with the analog output from a MZI signal as

---

schematically illustrated in Figure 4. Designed with OCT applications in mind, the AlazarTech series of digitizers are robust to the changing frequencies of the pixel clock signal, and signal dead time in between sweeps.

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For the cases where software resampling is preferred, the unused channel can be used to acquire a recalibration signal. The OCT data can then be resampled to be linear in  $k$ -space by using any of several algorithms (peak detection, zero crossings, etc.) on the recalibration signal. Both channels on the AlazarTech boards are digitized simultaneously, providing zero-shift between the OCT signal and the recalibration signal. For applications which require more than two inputs, multiple boards can be arranged in a Master/Slave configuration.

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AlazarTech has a wide-range of digitizers suitable for OCT applications. They range from 8-16 bit resolution cards with maximum sampling rates of 50MS/s to 1GS/s. For OCT applications, it is recommended to use at least 12-bit resolution digitizer cards. However, with advancement of high-speed SS OCT systems, there is a need for faster PC transfer rates. In this case, an 8-bit card (ATS9870 – 1GS/s, 8-bit) may be more preferred, which has been shown to only reduce the SNR of the OCT signal by  $\sim 0.6\text{dB}$  [7]. The wide range of AlazarTech waveform digitizers are complemented by an impressive library of examples and functions with the C/C++ and LabView Software Development Kits (SDK). OCT data acquisition and image processing can be easily integrated into the existing examples in the SDK.

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## References 参考书目

- [1] D. Huang, E. A. Swanson, C. P. Lin, J. S. Schuman, W. G. Stinson, W. Chang, M. R. Hee, T. Flotte, K. Gregory, C. A. Puliafito, and J. G. Fujimoto, "Optical Coherence Tomography," *Science* **254**, 1178-1181 (1991).
- [2] M. Wojtkowski, "High-speed optical coherence tomography: basics and applications High-speed optical coherence tomography: basics and applications," *Applied Optics*, 49 (16), D30-61, 2010.
- [3] J. Fujimoto, W. Drexler, *Optical Coherence Tomography: Technology and Applications*. Springer Berlin Heidelberg, 2008. [E-Book Series] Available: SpringerLink-Book.
- [4] B. E. Bouma, and G. J. Tearney, *Handbook of Optical Coherence Tomography*. Dekker, New York, NY, 2001.

- 
- [5] A. M. Rollins and J. A. Izatt, "Optimal interferometer designs for optical coherence tomography," *Optics Letters*. **24**, 1484-1486 (1999).
- [6] M. A. Choma, K. Hsu, J. A. Izatt, "Swept source optical coherence tomography using an all-fiber 1300-nm ring laser source," *Journal of Biomedical Optics* **10**(4), 044009 (2005).
- [7] B. D. Goldberg, B. J. Vakoc, W. Oh, M. J. Suter, S. Waxman, M. I. Freilich, B. E. Bouma, and G. J. Tearney, "Performance of reduced bit-depth acquisition for optical frequency domain imaging," *Opt. Express* **17**, 16957-16968 (2009).
- [1] D. Huang, E. A. Swanson, C. P. Lin, J. S. Schuman, W. G. Stinson, W. Chang, M. R. Hee, T. Flotte, K. Gregory, C. A. Puliafito, J. G. Fujimoto, "光学相干断层扫描", *科学*, 第 254 卷, 1178-1181 (1991).
- [2] M. Wojtkowski, "高速光学相干断层扫描: 基础和应用", *应用光学*, **49** (16), D30-61, 2010.
- [3] J. Fujimoto, W. Drexler, *光学相干断层扫描: 技术和应用*. 海德堡柏林: 施普林格出版社, 2008. [电子图书系列]可提供: SpringerLink 电子丛书.
- [4] B. E. Bouma, G. J. Tearney, *高速光学相干断层扫描手册*. 纽约州纽约市: 德克尔出版社, 2001.
- [5] A. M. Rollins, J. A. Izatt, "光学相干断层扫描的最佳干涉仪设计", *光学快报*, 第 24 卷, 1484-1486 (1999).
- [6] M. A. Choma, K. Hsu, J. A. Izatt, "使用全光纤 1300-nm 环形激光源的扫频光源光学相干断层扫描", *生物医学光学杂志*, 第 10(4)期, 044009 (2005).
- [7] B. D. Goldberg, B. J. Vakoc, W. Oh, M. J. Suter, S. Waxman, M. I. Freilich, B. E. Bouma, G. J. Tearney, "用于光频域成像的降低位深度采集的性能", *光学快讯*, 第 17 卷, 16957-16968 (2009).